

# JUVASONIC® GETTING STARTED

Use this guide to confidently deliver a complete procedure with precision, consistency, and proven results.



## Why Juvasonic®?

- » Enhances topical absorption
- » Delivers exosomes, growth factors, peptides, PRP + PDRN
- » Needle-free and pain-free,
- » No downtime
- » Efficient use of high-grade topicals

## Step-by-Step Instructions

1. Charge your Juvasonic® device.
2. Set up clean treatment area with all supplies.
3. Clean treatment area with alcohol and hypochlorous acid. Let dry.
4. Pre-condition skin/scalp with empty cartridge using vertical, horizontal, and circular motions for ~15-20 seconds per 2x2 area.
5. Withdraw serum with syringe. Dilute with saline, if needed.
6. Fill Juvasonic® 3 mL vial + place nanotip.
7. Insert cartridge. Select speed. Activate red/blue LED (double tap).
8. Apply serum using same motion technique, working in sections.
9. Finish with serum, optional mask, and sunscreen.

## Juvasonic Checklist

### Pre-Treatment

- Alcohol wipes
- Chlorhexidine
  - Gloves
  - 4x4 gauze
- Headband/towel

### Serum Delivery

- Juvasonic® (charged)
- 3 mL vial + nanotip
- 3 cc syringe
- 18G blunt tip needle
- Sterile saline
- Serum (Exosomes, PDRN, PRP, PRF, etc.)



### Post-Treatment

- Hyaluronic acid
- Soothing mask (optional)
- SPF 30+ sunscreen

### Contact Us

02 9528 4444  
info@enesi.com.au  
www.enesi.com.au  
www.juvasonic.com.au

